

Operating Instructions

Compact scale Puro® LargeCount



98628-000-63

Edition 1.0.2

6/25/2020

Foreword

Must be followed!

Any information in this document is subject to change without notice and does not represent a commitment on the part of Minebea Intec unless legally prescribed. This product should only be operated/installed by trained and qualified personnel. In correspondence concerning this product, the type, name, and release number/serial number as well as all license numbers relating to the product have to be cited.

Note

This document is partially protected by copyright. It may not be changed or copied, and it may not be used without purchasing or written permission from the copyright owner (Minebea Intec). The use of this product constitutes acceptance by you of the abovementioned provisions.

Table of contents

1	Introd	uction	4
1.1	Rea	d the manual	4
1.2	This	is what operating instructions look like	4
1.3	This	i is what lists look like	4
1.4	This	i is what menu items and softkeys look like	4
1.5	This	is what the safety instructions look like	4
2	Safety	instructions	6
2.1	Gen	eral information	6
2.2	lncc	oming goods inspection	6
2.3	B Befo	ore operational startup	6
	2.3.1	Danger of explosion	6
	2.3.2	IP protection	6
	2.3.3	Storage and transport conditions	7
2.4	Fail	ure and excessive stresses	7
3	Device	e installation	8
3.1	Med	hanical preparation	8
	3.1.1	Ambient conditions	8
	3.1.2	Installation location	8
	3.1.3	Unpacking	8
	3.1.4	Checking the equipment supplied	8
	3.1.5	Leveling the weighing platform	9
	3.1.6	Acclimatizing the device	9
3.2	2 Con	nection	9
	3.2.1	Electrical supply	9
	3.2.2	Connecting a printer	
	3.2.3	Connecting a scanner	
	3.2.4	Connecting a second weighing platform	
4	Device	e description	12
4.1	Disp	blay and operating elements	12
	4.1.1	Overview	12
	4.1.2	Display elements	13
	4.1.3	Operating elements	
5	Opera	ting	
5.1	Bas	ic functions	
	5.1.1	Switching on the device	
	5.1.2	Switching off the device	
	5.1.3	Adjusting the GEO setting	

	5.1.4	Display functions	17
	5.1.5	Increment d	17
	5.1.6	Scale with no load	17
	5.1.7	Enter values	17
	5.1.8	Initializations	17
	5.1.9	Applications activated during initial commissioning	17
	5.1.10	Connecting a second scale	17
5.2	Appl	ication programs	
	5.2.1	General information	
	5.2.2	Weighing application	
	5.2.3	Counting application	20
	5.2.4	Checking application	23
	5.2.5	Totalizing application	29
	5.2.6	Automatic tare	
	5.2.7	Automatic Printing	
	5.2.8	Product memory	
	5.2.9	Handling of the second scale	41
5.3	Men	u	41
	5.3.1	Menu display in displays	42
	5.3.2	Menu mode	44
	5.3.3	Processing entered values in the menu	45
	5.3.4	Menu navigation	46
5.4	Calib	ration, adjustment	54
	5.4.1	Calibration	54
	5.4.2	Linearization	57
	5.4.3	Geographic adjustment factor (GEO) – procedure	58
	5.4.4	GEO code table	59
5.5	Print		61
	5.5.1	Output format	62
	5.5.2	Product memory	63
	5.5.3	Menu settings: Print	63
5.6	PC o	utput	64
	5.6.1	SBI interface	64
6		nance/repairs/cleaning	
6.1		irs	
6.2		ning	
	6.2.1	Instructions for cleaning	
	6.2.2	Cleaning agents	66
7	Waste	disposal policy	
8		orrection	
8.1	Serv	ce information	69

9	Technical data	70
9.1	1 Specification	70
9.2	2 Accessories	71
	3 Dimensions	
	Appendix	
	Appendix	
10.1		75

1 Introduction

1.1 Read the manual

- Please read this manual carefully and completely before using the product.
- This manual is part of the product. Keep it in a safe and easily accessible location.

1.2 This is what operating instructions look like

1. - n. are placed before steps that must be done in sequence.

- is placed before a step.
 - ▷ describes the result of a step.

1.3 This is what lists look like

indicates an item in a list.

1.4 This is what menu items and softkeys look like

[] frame menu items and softkeys.

Example:

[Start]- [Applications]- [Excel]

1.5 This is what the safety instructions look like

Signal words indicate the severity of the danger involved when measures for preventing hazards are not followed.

▲ DANGER

Warning of personal injury

DANGER indicates death or severe, irreversible personal injury which will occur if the corresponding safety measures are not observed.

• Take the corresponding safety precautions.

△ WARNING

Warning of hazardous area and/or personal injury

WARNING indicates that death or severe, irreversible injury may occur if appropriate safety measures are not observed.

• Take the corresponding safety precautions.

Warning of personal injury.

CAUTION indicates that minor, reversible injury may occur if appropriate safety measures are not observed.

• Take the corresponding safety precautions.

NOTICE

Warning of damage to property and/or the environment.

NOTICE indicates that damage to property and/or the environment may occur if appropriate safety measures are not observed.

• Take the corresponding safety precautions.

Note:

User tips, useful information, and notes.

2 Safety instructions

2.1 General information

- The device may only be used as intended for weighing tasks.
- Observe the operating limits of the device.
- Do not apply loads that exceed the capacity of the scale.
- The voltage rating printed on the power supply (see type plate) must be the same as the local line voltage.
- Before connecting or disconnecting electronic peripheral devices, disconnect the device from the mains or from the data interface.
- Unplug the power cord from the mains supply before cleaning.
- Make sure that no liquid enters the device.
- The device may only be opened by authorized technicians.

2.2 Incoming goods inspection

Check the contents of the consignment for integrity. Check the contents visually to determine whether any damage has occurred during transport. If there are grounds for rejection of the goods, a claim must be filed with the carrier immediately. A Minebea Intec sales or service organization must also be notified. Visit our website http://www.puroscales.com or contact your dealer.

2.3 Before operational startup

NOTICE

Perform visual inspection.

Before operational startup as well as after storage or transport, inspect the product visually for signs of mechanical damage.

The product should not be put into operation if it displays signs of visible damage and/or is defective.

2.3.1 Danger of explosion

Do not use the device in hazardous areas.

2.3.2 IP protection

The model fulfills protection grade IP43.

2.3.3 Storage and transport conditions

NOTICE

Material damage is possible.

Unpacked devices may lose their precision due to strong vibrations; strong vibrations may impair the safety of the device.

Do not subject the device to extreme temperatures, moisture, shocks, and vibrations.

2.4 Failure and excessive stresses

If the device or the power cord display visible damage: Disconnect the power supply and secure the device to prevent it being used further.

Do not unnecessarily subject the device to extreme temperatures, corrosive chemical vapors, moisture, shocks, and vibrations.

Extreme electromagnetic influences can affect the display value. Once the disturbance has ceased, the product can be used again as intended.

3 Device installation

3.1 Mechanical preparation

3.1.1 Ambient conditions

- Only use within buildings.
- Operating temperature: -10°C to +40°C
- Storage temperature: -20°C to +50°C
- Relative humidity: 20% to 85%, non-condensing
- Altitude: up to 3,575 m

3.1.2 Installation location

- Place the device on a stable, flat surface.
- Position the device so that the power plug is freely accessible and the power cord does not present an obstacle or trip hazard.

Avoid unsuitable influences at the installation location:

- Extreme temperatures and excessive temperature fluctuations
- Heat due to proximity to heaters or due to direct sunlight
- Aggressive chemical vapors
- Extreme moisture
- Extreme vibrations

3.1.2.1 Shock resistance

NOTICE

Falling objects, side impacts, and shock loads may affect the performance and the accuracy of the scale and damage the platform.

Avoid shock loads!

3.1.3 Unpacking

- Unpack the device and check it for visible external damage.
 - ▷ If there is damage, follow the instructions in the chapter "Safety check".
- ► Keep the original packaging in case the device needs to be returned. Remove all cables before sending.

3.1.4 Checking the equipment supplied

- 1 scale
- 1 load plate
- 1 USB power supply with cable
- Safety instructions and QR code for access to the complete documentation


3.1.5 Leveling the weighing platform

To achieve reproducible weighing results at all times, the weighing platform must be set up to be precisely horizontal.

Therefore the weighing platform must be re-leveled every time it is moved to a different location.

Leveling the weighing platform

- Use the adjustable feet to align the weighing platform so that the air bubble of the level indicator is in the center of the circle.
- Check that all four of the adjustable feet are touching the surface.
 - ▷ The weight of the platform must be spread equally across the adjustable feet.
- Adjust the adjustable feet: Retract the adjustable feet (clockwise) in order to lift the scale. Extend the adjustable feet (counter-clockwise) in order to lower the scale.



3.1.6 Acclimatizing the device

If a cold device is brought into a warm environment, condensation may form.

• Keep the device disconnected from the mains and allow it to acclimatize at room temperature for approx. two hours.

3.2 Connection

3.2.1 Electrical supply

The scale is supplied using a power supply unit, unless a battery supply is required. Connect the USB-C male plug connector with the USB-C female plug connector on the underside of the device, then connect the power supply unit to the wall socket.

Note:

Do not use the USB-C connector cable for the PC communication. Instead, use a standard USB-C cable.

3.2.1.1 Battery power

The scale can be operated immediately with the power supply. In order to operate the scale with the battery, the battery should first be charged for 12 hours. If there is a power outage or if the power cord is disconnected, the scale switches into battery operation automatically. In the event of supply via a power supply, the battery is constantly charged meaning that the battery charging display (see Chapter 4.1.2) is continuously illuminated. The scale can be used during the charging process; the battery is protected against excess charging.

When the device is switched on, the battery status LED illuminates in red while the battery is charging, and it goes green when the battery is fully charged.

The battery must be charged in a dry environment. For a maximum operating time, the battery should be charged at room temperature.

During battery operation, the battery icon displays the battery's remaining charge status. The display switches off automatically when the batteries are empty.

lcon	Charge status
	0 to 10% remaining
	11 to 40% remaining
	41 to 70% remaining
	71 to 100% remaining

Note:

If the battery icon flashes rapidly, then there is around 30 minutes of working time left. When [lo.bat] is displayed, the scale switches off.

△ WARNING

Danger of explosion

If the rechargeable battery is replaced with a battery of the wrong type, or if it is not connected correctly, then there is a danger of explosion.

- The battery may only be replaced with the same type by an authorized Puro® service dealer.
- The battery must be disposed of according to the locally valid laws and regulations.

If the hardware does not recognize a connected rechargeable battery, the following applies:


Rechargeable battery	USB-C cable	LED display charging	Display icon
Battery not full	Connect	Red	No icon
Battery full	Connect	Green	No icon
No battery installed	Connect	Red	No icon
Battery not full	Disconnect	Switched off	Charge status of the battery
Battery full	Disconnect	Switched off	Battery full

3.2.2 Connecting a printer

A printer can be connected via the printer port on the underside of the device.

3.2.3 Connecting a scanner

Connect the scanner according to the following diagram. The scanner must be connected using the original cable YP-CYSSR1.




3.2.4 Connecting a second weighing platform

An analog weighing platform can be connected to the compact scale Puro® LargeCount and used as a second scale (weighing point).

- ► Turn the device on its side.
- Open the flap on the bottom of the device.
- Connect the weighing platform's 9-pin D-Sub male plug connector with the 9-pin D-Sub female plug connector. Pay attention to the pin assignment of the female plug connector and the color allocation of the cable.

Pin assignment of the 9-pin D-Sub female plug connector:



1	DGND	Ground
2	-	Not assigned
3	-	Not assigned
4	Ex+	Supply voltage+
5	Ex-	Supply voltage-
6	Se+	Sense+
7	Si+	Signal+
8	Si-	Signal-
9	Se-	Sense-

The A/D converter must then be configured and a two-point calibration carried out. For more information about the second weighing platform, see Chapters 5.1.10 and 5.3.4.9.


4 Device description

4.1 Display and operating elements


4.1.1 Overview

Display definitions:

- "Weight" (weight) display = small display on the left-hand side
- "PW" (sample weight) display = small display on the right-hand side
- "Count" (Counting) display = larger display in the middle




Control panel




No.	Description
1	Display elements, see Chapter <mark>4.1.2</mark> .
2	Operating elements, see Chapter 4.1.3.

4.1.2 Display elements




LCD display



ltem	Description	ltem	Description
1	Busy (process running)	2	Plus sign
3	Minus sign	4	1/4 d range around zero
5	Counting application active	6	Weighing application active
7	Checking application active	8	Totalizing application active
9	Icon for automatic tare or automatic reference optimization	10	Average sample weight too low
11	Sample weight too low	12	Two scales are active
13	Number of the active scale	14	Data transmission active
15	Weight unit	16	Battery charging
17	Warning icon	18	Printer icon
19	Net, preset, tare value, gross	20	Item (value in items)

LED displays



The colored LED displays (1) on the right-hand side of the control panel for counting are used in the Checking application (see Chapter 5.2.4.2.2) and illuminate according to the following rules:



(Red) weight value > upper tolerance limit (Green) weight value $\leq \geq$ lies within the tolerance limits (Yellow) weight value < lower tolerance limit

4.1.3 **Operating elements**

Key	Primary function (Brief press) < 1 second	Secondary function (Extended press) holding > 2 seconds
On/Zero Off	On/Zero Switch on the scale (if the scale is switched off), zero (if the scale is switched on)	Off Switch off the scale
Print Unit	Print Send the current value to the selected COM ports if the "Out" option is specified for auto- matic printing.	Unit Change weight unit
Target	Target Specify tolerance limits if the Check Weighing or Check Counting application is active	Select Check Weighing applicati- on
M+ Menu	M+ Write into the totalizing memory or exit the entry process.	Menu Call up the user menu
Mem	Mem Save or load a product	Start product defini- tion
Sample	Sample Adopt quantity of the reference weights and calculate the sample weight	
APW	Average sample weight Adopt sample weight	

Key	Primary function (Brief press) < 1 second	Secondary function (Extended press) holding > 2 seconds
)→ T ←	Tare Specify or delete a tare value	Delete totalizing me- mory
9	09 Enter characters	
* 1 2	*12	Switch scale
C +/-	C +/- Delete the last character entered or exit the dis- play of the totalized values	Switch the algebraic sign

Keys for the menu navigation

Key	Menu function (Brief press) < 1 second
On/Zero Off	Yes Adopt the current setting in the display
Print Unit	No Discard the current setting in the display and switch to the next availa- ble setting Switch to the next menu or item/display the next value
Target	Back Switch to the previous menu items/display the previous value
M+ Menu	Exit Exit the user menu Cancel the ongoing calibration


5 **Operating**

5.1 Basic functions

5.1.1 Switching on the device

- ► Briefly press the On/Zero Off key.
 - > All elements of **all** displays are displayed for 2 seconds.

All LEDs of the Checking application illuminate for 2 seconds.



Then the software version number (here [SR 2.5]) is displayed for 2 seconds in the "Weight" display; at the same time, the selected GEO area is displayed in the "Piece Weight" display, e.g.: [GEO 12]. Nothing is displayed in the "Count" display.

Weight

Piece Weight



The combination of applications that was most recently active (selected) before switching off is started. If the applications had not yet been initialized, they start with the following parameters:


The default unit when starting for the first time is [kg].

If a second scale is defined in the menu, the device always starts with WP1 (weighing platform 1).

5.1.2 Switching off the device

- Press and hold the On/Zero Off key.
 - ▷ [-OFF-] is displayed for approx. 2 seconds in the "Count" display.

Count



The other displays are dark. Then the device switches off and all displays are dark.

This key is active in **ALL** application and menu statuses.

5.1.3 Adjusting the GEO setting

Adjust the GEO setting according to the location in order to guarantee accurate weighing results. See Chapter **5.4.3**.

5.1.4 Display functions

In the normal weighing mode, the displays have the following functions:

- The actual weight value is displayed with the selected unit in the "Weight" display.
- The actual sample weight in grams or ounces is displayed in the "PW" display.
- The calculated number of parts is displayed with the unit [pcs] (parts) in the "Count" display.

The icon for "Charging" 🖹 is only displayed in the "Weight" display if a rechargeable battery is connected.

The icon for "active scale" $\overline{\Delta\Delta}$ 12 is only displayed in the "Weight" display. If only one weighing platform is connected, the icon for "active scale" is dark.

5.1.5 Increment d

"d" stands for the lowest weight value that can be displayed.

Example d = $0.02 \text{ g} \rightarrow 2 \text{ d} = 0.04 \text{ g} \rightarrow 3 \text{ d} = 0.06 \text{ g}$

5.1.6 Scale with no load

In general there is no load on a scale if the weight on the weighing pan is under 2 d (< 2 d).

5.1.7 Enter values

Values can be entered via the keypad; the display takes place in the "Count" display. At the start of an entry, the display is dark and on the right-hand side a flashing

underscore (cursor) marks the last position. The entry can be corrected using the key, which is used to delete the last digit of the entry in each case.

Alternatively, the decimal point $\begin{bmatrix} * \\ 1 \\ 2 \end{bmatrix}$ can also be entered as the first character. The software then automatically adds the zero in front of the decimal point.

5.1.8 Initializations

During an initialization it is not possible to activate the menu.

5.1.9 Applications activated during initial commissioning

The Weighing and Totalizing applications are automatically activated during the initial commissioning of the device. Counting is active, but not initialized. The Checking application is switched off.

5.1.10 Connecting a second scale

A second scale can be connected in order to achieve a counting system that consists of a reference scale and a scale for larger amounts.

Each scale has a separate serial number with 13 characters.

Number concept:

- Serial number for scale 1: W1 38457989
- Serial number for scale 2: W2 38457989

The following is printed:

```
SERNO: W1 38457989
SERNO: W2 38457989
```

The following is displayed under [INFO] in the menu: [W1 38457989] and [W2 38457989]

5.2 Application programs

5.2.1 General information

In principle the main Counting application is always activated.

The following applications can be activated in parallel to the Counting application:

- Checking (two applications: Check Weighing and Check Counting)
- Totalizing (for totalizing parts or weights)
- Automatic tare
- Automatic printout

The product memory is available for 30 products.

5.2.2 Weighing application

If the Counting application has not been initialized (sample weight = 0), the application

```
icon for weighing 4 is displayed in the "Weight" display.
```

Weight

5.2.2.1 Stability

The weight value and the application icons are displayed in the "Weight" display.

If a weight value is stable, the unit symbol is displayed in the "Weight" display and in the "Count" display:

stable:	[2.342 kg]	and [47 pcs]
ot stable:	[2.342]	and [47]

and printed:

```
stable:
not stable:
```

2.342 kg, 47 pcs 2.342 , 47

Stable weight value in the "Weight" display:

Weight

•	

5.2.2.2 Set tare

- ► To tare, briefly press the $\rightarrow T \leftarrow$ key when there is a load on the weighing pan.
 - ▷ The icon [NET] is displayed in the "Weight" display.

During the tare process, the busy icon is displayed in the "Weight" display with no weight value; the "Count" display is dark.

Tared value in the "Weight" display:



5.2.2.3 Preset tare value

There is only one tare memory in the device, which contains either a measured

weight value (press $\rightarrow T \leftarrow$ key) or an entered value (preset tare value).

- 1. Enter a weight value such as [0.010] via the numeric keypad or using a scanner.
 - ▷ This value is displayed in the "Count" display.

Count



2. Briefly press the $\rightarrow T \leftarrow$ key.

▷ The entered value is specified as the preset tare value with the corresponding unit; the value currently in the tare memory is overwritten.

The weight value in the "Weight" display is calculated in advance using this preset tare value.

The actual number of parts on the weighing pan is recalculated and displayed in the "Count" display.

The preset tare value can be deleted as follows:

3. Press the $\rightarrow T \leftarrow$ key if there is no load on the weighing pan.

 \triangleright Delete the tare memory.

Specify a new tare value as follows:

Press the $\rightarrow T \leftarrow$ key, if there is a load on the weighing pan, in order to overwrite the tare memory with the weight on the weighing pan.

Or press the Or/Zero of the tare memory. The preset tare value is deleted after the device is switched off.

5.2.2.4 Weight unit key

The display of a weight value can be switched between various weight units.

Pressing and holding the
Print Unit key changes the displayed unit.

Possible units are: Gram [g], Kilogram [kg], Pound [lb], Ounce [oz], Pound-ounce [l:o].

Only units that have been activated in the menu can be changed.

While the key is held down, the unit changes and after around 2 seconds the next unit is displayed, and so on, until the key is released.

This is how the switched units are displayed in the "Weight" display and in the "PW" display:

Units in the "Weight" display	Units in the "PW" display
Kilogram	Gram
Gram	Gram
Pound	Ounce
Ounce	Ounce
Pound: Ounce	Ounce

5.2.2.5 **Printouts**

The elements to be printed are configured in the menu.

Printout		Description
5.003 g	N	Positive net weight value
- 0.003 g	N	Negative net weight value
2.003 g	G	Positive measured gross weight value
2.003 g	G#	Positive calculated gross weight value
1.003 g	Т	Tare weight value (measured value)
0.010 g	ΡT	Preset tare value (entered value)
- 0.010 !	G	Stable gross weight below zero

5.2.3 **Counting application**

The scale starts with "Counting". If a sample weight had been specified before switching off the device, the device starts with this sample weight.

The Counting application is always active. However, there is the option for this application not to be initialized.

If no sample weight has been initialized, both displays show [0.] and the application icon শ

(weighing) is displayed in the "Weight" display.

The device starts with the Totalizing application activated. Weight





Counting is initialized

5.2.3.1 Initializing the Counting application

- Initialization with a known reference quantity (see Chapter 5.2.3.1.1)
- Initialization with a known sample weight (see Chapter 5.2.3.1.2)

The initialization of the Counting application is only possible on scale 1. The initialization is carried out with the internal scale resolution.

5.2.3.1.1 Initialization with known reference quantity

This initialization is only carried out when the weight values are stable.


- Initialization (see Chapter 5.2.3.1.1.1)
- Initialization of a negative sample (see Chapter 5.2.3.1.1.2)

These initializations are possible even if a sample weight has already been calculated.

The entry can be corrected using the $\begin{bmatrix} c \\ +/- \end{bmatrix}$ key, which is used to delete the last digit of the entry in each case.

If there is no load on the weighing pan (< 2 d), the error message [REF.ERR] is displayed for around 2 seconds in the "PW" display.

Piece Weight



5.2.3.1.1.1 Initializing

- 1. Place some parts on the weighing pan.
- 2. Enter the known number of parts on the weighing pan via the keypad (only whole number values) or using a scanner.
 - ▷ The entered value is displayed in the "Count" display.
- 3. Briefly press the key.
 - ▷ The newly calculated sample weight is displayed in the "PW" display, and the actually calculated number of parts is displayed in the "Count" display.

5.2.3.1.1.2 Initialization via the removal of parts

If the removal of parts from a container is required:

- 1. Place a full container on the scale.
- 2. Briefly press the $\rightarrow T \leftarrow$ key (tare container).
- 3. Remove some parts from the container.
 - ▷ Now a negative weight value is displayed in the "Weight" display.

Proceed with step 2 in Chapter 5.2.3.1.1.1.

5.2.3.1.2 Initialization according to known sample weight

- 1. Enter the desired sample weight via the keypad or using a scanner.
 - ▷ The entered value is displayed in the "Count" display.
- 2. Briefly press the key.
 - ▷ The newly adopted sample weight is displayed in the "PW" display, and the actually calculated number of parts is displayed in the "Count" display.

5.2.3.2 Performing Counting (without check weighing)

Automatic reference optimization (see Chapter 5.2.3.2.1).

If a sample weight is initialized, it is displayed in the "PW" display.

- Place a weight on the weighing pan.
 - ▷ The actual weight value is displayed in the "Weight" display.

The calculated number of parts is displayed in the "Count" display.

The application icon $\delta \delta$ (Counting) is displayed in the "Weight" display.

- Briefly pressing the Unit Unit key will generate a printout.
- Pressing and holding the Print Unit Unit key changes the unit in the "Weight" display (in certain cases also in the "PW" display [see Chapter 5.2.4.1]).
- Briefly pressing the key or the sample key or the keypad starts an initialization as described above.
- Press and hold or briefly press the key (see Chapter 5.2.4)
- Press the Mem key (see Chapter 5.2.8)
- Pressing the +/key will delete the sample weight.

5.2.3.2.1 Automatic reference optimization

If in the menu [OP.FUNC]- [A.OPT]- [ON] is set, the automatic reference optimization is activated. The icon is displayed in the "PW" display.

The main feature is the automatic more precise calculation of the sample weight.

The sample weight is recalculated during the optimization if the following conditions are all met:

- Scale is at a standstill
- No prefix change
- The current quantity exceeds the original reference quantity by at least two. The new quantity may not exceed twice the old reference quantity (this limitation does not apply for the first optimization if the sample weight was entered using a bar code scanner or the keyboard).
- The internally calculated quantity (e.g. 17.24) must be less than ± 0.3 parts different from the whole number (in the example: 17).
- AUTO - If an optimization has been carried out, the icon OPT is displayed in the "PW" display, and the newly calculated sample weight is also displayed in the "PW" display.

The automatic optimization is not carried out on scale 2.

5.2.3.3 Printouts (without checking)

Normal printout:

Printout	1		Description
	441 pcs	QNT	Positive value
	- 41 pcs	QNT	Negative value
MODE:	COUNT		Activated application
WREF	4.15431 o	Z	Sample weight, same value as that displayed in the
			"PW" display

5.2.4 Checking application

For the Checking application, there are indicator LEDs that show the current range.

5.2.4.1 Activating the checking application

While the Counting application is being performed, the checking applications can be activated via the $\begin{bmatrix} Target \\ Target \end{bmatrix}$ key.

Press and hold the key.

▷ [CHECK] is shown in the "Weight" display.

Weight

Piece Weight



The activated checking mode is displayed for approx. 2 seconds in the "PW" display. The "Count" display is dark.

After this waiting time the next mode is displayed and so on. If the key is released, the mode displayed in the "PW" display is selected.

The following selection is possible:

[OFF]	Checking is switched off.
[WEIGHT]	Check Weighing is activated
[COUNT]	Check Counting is activated

If the Check Weighing application [WEIGHT] is activated and the Counting application has not been initialized, the following application icons will be $\Delta T \Delta + I$

displayed in the "Weight" display: 4 +/-.

If the Counting application [COUNT] has been initialized, the Counting $\Delta T \Delta + I$

If the Check Counting application [COUNT] is activated and the Counting application has been initialized, the following application icons will be displayed in the "Weight" display: $\circ \circ \circ$ +/-.

5.2.4.2 Checking weight values

- Initializing Checking application (see Chapter 5.2.4.2.1)
- Execution mode (see Chapter 5.2.4.2.2)
- Print outs (see Chapter 5.2.4.2.3)

5.2.4.2.1 Initializing checking

Pressing and holding the $\begin{bmatrix} C \\ +/- \end{bmatrix}$ key will generate a minus algebraic sign in the entry mode in the "Count" display.

Pressing and holding the $\frac{M^+}{Menu}$ (Exit) key will stop the initialization process immediately, without saving a new entry or limit.

Briefly pressing the Print Unit (No) key or Print Unit (No) key or Print Unit will delete the flashing limit value; on the right-hand side, a flashing cursor indicates that a new entry can be made. Initialization is possible on both scales.

- 1. After selecting a checking application, briefly press the key.
 - ▷ [SET.LOW] is shown in the "Weight" display.

Veight			
5	E	T .L	

The previous lower limit is also displayed in the "PW" display with the unit that

was activated before pressing the <u>key</u>.

This weight value is also displayed flashing in the "Count" display. The yellow LED illuminates.

- 2. Use the keypad or a scanner to enter a numerical weight value, which is then displayed in the "Count" display.
 - ▷ The "Count" display shows the same unit as in the "Weight" display.

The entered value does not flash. Only the actual cursor position (which is designated by an underscore) flashes. The entry starts on the right-hand side.

- 3. Confirm the entry by briefly pressing the Original (Yes) key.
 - Briefly pressing the (Back) key will discard the entry; the limit is not changed.

[SET.HI] is shown in the "Weight" display.

Weight

5E T.H I	56	Ţ. '.	
----------	----	--------------	--

The previous upper limit is also displayed in the "PW" display with the unit that Target

was activated before pressing the key.

This weight value is also displayed flashing in the "Count" display. The red LED illuminates.

- 4. Use the keypad or a scanner to enter a numerical weight value, which is then displayed in the "Count" display.
 - ▷ The entered value does not flash. Only the actual cursor position (which is designated by an underscore) flashes.
- 5. Confirm the entry by briefly pressing the Orff (Yes) key or discard by pressing the key.
 - ▷ If the limit conditions are correct, the checking application is initialized.

Count



If there is an error present, because, e.g., upper limit < lower limit, [LIM.ERR] is displayed for roughly 2 seconds in the "Count" display and the display returns to step 1.

The normal Counting application is activated again. In addition, the test LEDs are activated when there is a load on the weighing pan.

When checking additive weights, both limits must be a positive value.

When checking removed weights, both limits must be a negative value.

When using "Check Against Zero", the upper limit is a positive value and the lower limit is a negative value.

When checking a precise weight, the lower and upper limits must be the same.

5.2.4.2.2 Execution mode

Limit value LEDs:

positive weight	< lower limit	Yellow LED illuminates
positive weight	\geq lower limit and \leq upper limit	Green LED illuminates
positive weight	> upper limit	Red LED illuminates
negative weight	> lower limit	Yellow LED illuminates
negative weight	\leq lower limit and \geq upper limit	Green LED illuminates
negative weight	< upper limit	Red LED illuminates

When using "Check Against Zero", the reference weight is placed on the weighing pan and then the scale is tared. The reference weight is removed; then the weight to be checked is placed on the weighing pan.

5.2.4.2.3 Printouts

Counting is uninitialized:

Printout	Description
115 g OVER	Positive net value > upper limit
- 115 g OVER	Negative net value < negative upper limit
99 g ACCEPT	Positive net value in the target range
75 g UNDER	Positive net value < lower limit
MODE: CHECKWEIGH	Activated application
UNDER LIMIT 81 g	Lower limit
OVER LIMIT 100 g	Upper limit

Printout	Description
115 g OVER	Positive net value > upper limit
- 115 g OVER	Negative net value < negative upper limit
99 g ACCEPT	Positive net value in the target range
75 g under	Positive net value < lower limit
115 pcs	Net quantity
MODE: CHECKWEIGH	Activated application
UNDER LIMIT 81 g	Lower limit
OVER LIMIT 100 g	Upper limit
WREF 0.35423 g	Sample weight

Counting is initialized:

5.2.4.3 Checking number of parts

- Initializing Check Counting (see Chapter 5.2.4.3.1)
- Execution mode (see Chapter 5.2.4.3.2)
- Print outs (see Chapter 5.2.4.3.3)

5.2.4.3.1 Initializing Check Counting

Pressing and holding the $\begin{bmatrix} C \\ +/- \end{bmatrix}$ key will generate a minus algebraic sign in the entry mode in the "Count" display.

Pressing and holding the $\frac{M^+}{Menu}$ (Exit) key will stop the initialization process immediately, without saving a new entry or limit.

Briefly pressing the $Print \\ Unit \\ Unit \\ (No)$ key or $c \\ +/-$ will delete the flashing limit value; on the right-hand side, a flashing cursor indicates that a new entry can be made. Initialization is possible on both scales.

- 1. After selecting a checking application, briefly press the key.
 - ▷ [SET.LOW] is shown in the "Weight" display.

Weight

- L	E	T .L	

The previous lower limit is displayed in the "PW" display.


This number of parts is also displayed flashing in the "Count" display. The yellow LED illuminates.

Target

- 2. Use the keypad or a scanner to enter a number of parts, which is then displayed in the "Count" display.
 - ▷ The entered value does not flash. Only the actual cursor position (which is designated by an underscore) flashes. The entry starts on the right-hand side.
- 3. Confirm the entry by briefly pressing the On/Zero Off (Yes) key.

▷ Briefly pressing the (Back) key will discard the entry \rightarrow the limit is not changed.

[SET.HI] is shown in the "Weight" display.



The previous upper limit is also displayed in the "PW" display.

This number of parts is also displayed flashing in the "Count" display. The red LED illuminates.

- 4. Use the keypad or a scanner to enter a number of parts, which is then displayed in the "Count" display.
 - ▷ The entered value does not flash. Only the actual cursor position (which is designated by an underscore) flashes.
- 5. Confirm the entry by briefly pressing the On/Zero Off (Yes) key or discard by pressing the key.
 - ▷ If the limit conditions are correct, the checking application is initialized.

Count



If there is an error present, because, e.g., upper limit < lower limit, [LIM.ERR] is displayed for roughly 2 seconds in the "Count" display and the display returns to step 1.

The normal Counting application is activated again. In addition, the test LEDs are activated when there is a load on the weighing pan.

When checking additive weights, both limits must be a positive value.

When checking removed weights, both limits must be a negative value.

When using "Check Against Zero", the upper limit is a positive value and the lower limit is a negative value.

When checking a precise weight, the lower and upper limits must be the same.

5.2.4.3.2 Execution mode

- Only the number of parts is checked.

Limit value LEDs:

Positive number of parts	< lower limit	Yellow LED illuminates
Positive number of parts	\geq lower limit and \leq upper limit	Green LED illuminates
Positive number of parts	> upper limit	Red LED illuminates
Negative number of parts	> lower limit	Yellow LED illuminates

Negative number of parts	\leq lower limit and \geq upper limit	Green LED illuminates
Negative number of parts	< upper limit	Red LED illuminates

When using "Check Against Zero", the reference quantity to be checked is placed on the weighing pan and then the scale is tared. The reference quantities are removed; then the number of parts to be checked is placed on the weighing pan.

The warning symbol Δ appears in the "Count" display if the unit [pcs] is displayed.

5.2.4.3.3 Printouts

Counting is uninitialized:

Printout	Description
115 g N	Normal net value as result
MODE: CHECKWEIGH	Activated application
UNDER LIMIT -20 pcs	Lower limit
OVER LIMIT -60 pcs	Upper limit

Counting is initialized:

Printout	Description
115 pcs OVER	Positive net number > upper limit
- 115 pcs OVER	Negative net number < negative upper limit
99 pcs ACCEPT	Positive net number in the target range
75 pcs UNDER	Positive net number < lower limit
MODE: CHECKCOUNT	Activated application
UNDER LIMIT -20 pcs	Lower limit
OVER LIMIT -60 pcs	Upper limit
WREF 0.35423 oz	Sample weight, as displayed

5.2.4.4 Menu settings

For the Checking application there are other functions that can be activated via the menu:

OP.FUNC	
– A.TARE	Automatic tare
- OFF	deactivated* (default setting)
– ON	1st stable weight is tared
– ON-ACC	Stable loads within the acceptance limits are tared (in all Checking applications)
– BEEP.SI	Signal (in the check weighing mode)
- OFF	deactivated* (default setting)
– ACCEPT	Signal when the weight is within the tolerance range
– UNDER	Signal when the weight is below the lower limit
– OVER	Signal when the weight is above the upper limit
– UNDOVR	Signal when the weight is outside of the tolerance range

5.2.5 Totalizing application

The application icon looks like this: Σ . It is only displayed in the "Weight" display. If the Counting application has not been initialized, only weights are totalized. If the Counting application has already been initialized, parts and weights are totalized. Only net values can be totalized.

5.2.5.1 Menu selection

This application can be selected in the menu:

OP.FUNC	
– TOT.SET	Totalizing setting
– OFF	deactivated* (default setting)
– AUTO	Automatic totalizing
— MAN	Manual totalizing


5.2.5.2 Execution modes

- Manual mode (see Chapter 5.2.5.2.1)
- Automatic mode (see Chapter 5.2.5.2.2)
- Print outs (see Chapter 5.2.5.2.3)
- The statistics information can be displayed by briefly pressing the display.
 M+
 Menu
 key in the

For this purpose there must be no load on the weighing pan (weight < 2 d).

The following parameters are displayed for three seconds:

Weight



Totalized weight with unit

Piece Weight

Number of items in the totalizing memory

Count



If counting is active, the total parts are

displayed. If only weighing is active, the display is dark.

Once the first three seconds have elapsed, the next statistics parameters are also displayed for three seconds:

Weight



Piece Weight	
	Min. weight
Piece Weight	······································
	Min. number of parts
Count	J
	Max. weight
Count	
pcs	
active.	Maximum number of parts if counting is
The $\begin{bmatrix} \mathbf{C} \\ +/- \end{bmatrix}$ key is used to exit this status p status.	prematurely and reactivate the previous
This status is exited automatically when the previous status takes place.	the waiting time has expired. The switch to
Statistics information can be deleted by	pressing and holding the $\rightarrow T \leftarrow$ key.
For this purpose there must be no load o pan < 2 d.	on the weighing pan: Weight on the weighing
Weight	
[CLR.TOT] is now displayed.	
If the statistics information needs to be	deleted, briefly press the On/Zero Off (Yes) key.
If the statistics information does not neekey.	

The initialization of Counting can be canceled by briefly pressing the key or by initializing Counting again; in this case the totalizing parameters are deleted.

After writing into the totalizing memory, the "Weight" display shows the total weight (= total value following adoption), the "PW" display shows the number of items in the totalizing memory, and the "Count" display shows the totalized number of parts. This display lasts for around 2 seconds.

If a weight or a number of parts has been adopted into the totalizing memory, the total symbol flashes until the weight is removed from the weighing pan and the weight is < 2 d. A new totalizing process can only be started if the load has previously been removed from the weighing pan.

If two weighing platforms are connected, the displayed weights are totalized accurately according to the display. Example: The first element in the totalizing memory is taken from WP1 and the second is taken from WP2:

1.003 kg	// WP1 with accuracy 1 g
= 5.15 kg	// WP2 with accuracy 50 g
= 6.153 kg	// total = totalizing memory

It is also possible to totalize negative weight values (deduction mode). Place weights on

the weighing pan, press the $\rightarrow T \leftarrow$ (Tare) key. Remove the first weight. Add to the

totalizing memory. Press the $\rightarrow T \leftarrow$ (Tare) key again. Then remove the second weight. Add to the totalizing memory.

5.2.5.2.1 Manual mode

Stable weights ≥ 2 d can be totalized by briefly pressing the \underbrace{Menu}_{Menu} key. This also applies for a number of parts if there is stability.

The next weight/the next quantity can be totalized if the load had previously been removed from the scale.

5.2.5.2.2 Automatic mode

Stable weights ≥ 2 d are automatically totalized if the load had previously been removed from the scale. This also applies for a quantity if there is stability.

5.2.5.2.3 Printouts

Min, Max is only printed if "All" is selected as print content in the menu [PRINT]-[CONTNT]- [TOTAL].


Weights as parameters (= weight total printout):

Printout		Description
N: 2		Number of items, here: 2
TOTAL:	1.955 g	Totalized value
MAX:	1.485 g	Maximum value
MIN:	0.470 g	Minimum value

Parts as parameters (= parts and weight total printout)

Printout		Description
N: 25		Number of items, here: 25
TOTAL:	148 g	Totalized value of the weights
	124 pcs	Totalized value of the parts
MAX:	20 g	Maximum value as weight
	10 pcs	Maximum value as number of parts
MIN:	4 g	Minimum value as weight
	2 pcs	Minimum value as quantity

5.2.6 Automatic tare



5.2.6.1 Menu selection

OP.F	UNC
– A. 1	FARE
	— OFF
	— ON
	– ON-ACC

Automatic tare deactivated* (default setting) first stable weight ≥ 2 d is tared Test: Each stable weight within the acceptance values is tared. This means that the automatic tare is not just carried out once!

5.2.6.2 Execution mode

Standard performance:

- The first weight that is placed on the weighing pan is tared if $\ge 2 d$.
- If there is no load on the scale (< 2 d), the tare memory is deleted.

Performance during "checking":

- First option [ON]:

The first weight ($\ge 2 d$) that is placed on the weighing pan is tared

Second option [ON-ACC]:

In the event of load ≥ 2 d, only stable values within the tolerance range are tared. Here, each weight in the accepted limit range is tared.

This function is carried out before the automatic printing and also before the automatic totalizing. The automatic tare has a higher priority than printing and totalizing. The first applied weight is tared, the second is adopted into the totalizing memory and automatically printed out if "Automatic Printing" and "Automatic Totalizing" are selected in the menu.

"Automatic tare" is not carried out if the menu is activated or an initialization is being carried out.

"Automatic tare" is deactivated if a preset tare value has been entered or a product with preset tare value has been activated.

5.2.7 Automatic Printing

Printing is carried out automatically **once** if the net weight on the weighing pan is ≥ 2 d. If the weight is removed from the weighing pan and the net weight is < d, the next net weight ≥ 2 d generates printing again.

5.2.7.1 Menu selection

Specify conditions under the menu item [PRINT] (printing):

PRINT	
– A.PRINT	automatic printing
- OFF	deactivated* (default setting)

— ON.STAB	automatic printing at standstill once, if weight > 1 d
— INTER	automatic printing in defined second intervals without standstill
- 25	Interval can be entered in seconds in the range from 1 to 3,600
– CONT	automatic printing for each weighing cycle without standstill
— ACCEPT	automatic printing once at standstill within the control limits

5.2.7.2 **Execution mode**

- Execution mode, general (see Chapter 5.2.7.2.1)
- Execution mode: ON STABLE (see Chapter) -
- Execution mode: INTER (see Chapter 5.2.7.2.2)
- Execution mode: CONTNT (see Chapter 5.2.7.2.3)
- Execution mode: ACCEPT (see Chapter 5.2.7.2.4) -

5.2.7.2.1 Execution mode, general


"Automatic Printing" normally has second priority after "Automatic tare". However, if "Automatic Totalizing" is also active then this is carried out before the printing. Each selected element in [PRINT]- [CONTNT] is printed, apart from the total parameters if the Totalizing application is not activated in the menu.

5.2.7.2.2 Execution mode: INTER

An interval is selected in the menu. The divisions can be selected in seconds from 1 to 3.600 seconds.

Example: If 5 seconds is selected, then a printout is created automatically every 5 seconds. The printout is carried out irrespective of whether or not the weight value is stable.

Weight



In this status, the number of seconds that can be entered via the keypad or using a scanner is displayed in the "Weight" [A.PRINT] display, in the "PW" [INTER] display, and in the "Count" display.

5.2.7.2.3 Execution mode: CONTNT

The elements are printed as quickly as possible with and without standstill.

5.2.7.2.4 Execution mode: ACCEPT

Automatic Printing is only carried out when the "Checking" application is activated AND the weight applied or the weight removed is within the control limits.

This automatic printing is generated **after the** totalizing (automatic mode), if "Totalizing" is activated.

5.2.8 **Product memory**

The device has a product memory for a maximum of 30 products.

Each product contains:

- Product ID (memory item number)
- Product name = 12 ASCII characters long. Alphanumeric characters are possible when using a scanner!
- Preset tare value
- Sample weight
- "Checking" lower limit = a weight value or a number of parts.
- "Checking" upper limit = a weight value or a number of parts.

If a checking limit is not equal to zero, the "checking" application is activated automatically. If these limits are zero, "checking" is automatically deactivated.

Example: "Checking" is activated. When a product is loaded that only contains counting parameters, "checking" is set to [OFF] (deactivated).

Mem key is pressed briefly and the user wants to exit this status (saving or If the M+ Menu (Exit) key.

loading), briefly press the

If two scales are defined in the menu, the product can still be loaded even if the scales have different accuracies. This is the customer's responsibility.

5.2.8.1 Saving product data

Requirement for saving: The application has already been initialized.

- Mem 1. Briefly press the key.
 - The following is displayed in the displays: \triangleright



The memory number (ID) can be changed using the keyboard or the scanner.

On/Zero Mem Off 2. Briefly press the (Yes) key in order to save the product parameters.

The first product memory has the number 1.

If the entry of the memory number contains an error, e.g. if the memory number is too big, then [LIM.ERR] is displayed in the "Count" display for approx. 2 seconds.

Count



During the saving process, the busy icon \bigotimes is activated in the Weight display; the other displays are dark.

The memory number can be freely selected. It does not have to be the next free one.

If the selected memory is already occupied, this memory is overwritten.

5.2.8.2 Loading product parameters

- 1. Enter a product memory number via the keypad or using a scanner.
 - ▷ The number is displayed in the "Count" display.
- 2. Briefly press the key.
 - ▷ The following is displayed in the displays:



- 3. Press the _____ or _____
 - ▷ The product parameters of the selected product are loaded.

(Yes) key.

If the memory is empty, [NO.DATA] is displayed for approx. 1 second in the "Weight" display.

Weight

Then the device switches back into the status from before entering the product number.

The product ID in the "Count" display can be changed by briefly pressing the $\underbrace{Print \\ Unit}$ (No)

and L_____ (Back) keys.

Pressing the Unit (No) key will incrementally increase the number.

Pressing the (Back) key will incrementally decrease the number.

Example: Memory items 1, 2, 5, 7 are occupied. If the displayed number is 5, pressing the

 $\frac{\frac{Print}{Unit}}{(No) \text{ key will display}} \rightarrow 7. \text{ Pressing} \qquad (Back) \text{ will display} \rightarrow 5.$

5.2.8.3 Changing/defining product parameters

A current product can be changed or a new product can be defined.

Briefly pressing the (Exit) key in any status of this entry process will reset the device into the status before starting this entry process.

- All entries can be made via the keyboard or using a scanner.
- The unit of the parameters is the one that was used before activating this change process.
- 1. Enter the desired product ID via the keypad or using a scanner.
 - ▷ The number is displayed in the "Count" display.
- 2. Briefly press the key.
 - ▷ The following is displayed in the displays:



- 3. Press and hold the key.
 - [CHANGE] is displayed in the "Weight" display for as long as the key is held down. Weight



If the key is released, [PROD.NA] is displayed in the "Weight" display for 2 seconds.

Piece Weight

|--|

After these 2 seconds, the saved product name is displayed flashing in the "Weight" display and in the "PW" display. The entered product ID is displayed in the "Count" display (no change).

If a product name is entered while [PROD.NA] is being displayed, the entry process starts directly with a flashing underscore (cursor) of the last digit in the "PW" display.

The product name, which may be up to 12 characters long, can be entered in the "PW" display and in the "Weight" display. The characters that can be used are numbers, capital letters, "-", "/", ".", and spaces.

Example: The product name is 123456789012.


Then [123456] is displayed in the "Weight" display and [789012] is displayed in the "PW" display.

- 4. Briefly press the $\begin{bmatrix} c \\ +/- \end{bmatrix}$ key while the product name is flashing.
 - The entire product name is deleted in the displays, and a flashing underscore is displayed at the last position for a new entry in the "PW" display.
- 5. Briefly press the $\begin{bmatrix} c \\ +/- \end{bmatrix}$ key while the cursor is flashing.
 - ▷ The last digit before the flashing cursor is deleted.

Pressing a key on the keypad will delete the old name and display the printed number.

- 6. Briefly press the On/Zero Off (Yes) key.
 - ▷ The entered sequence of characters is adopted. Continue with point ①.
- ① Display up

Weight



Piece Weight

Flashing weight value of the preset tare

value



Adopt the value by briefly pressing the On/Zero Off Off (Yes) key. Continue with point (2).


Press the (Back) key to return to entering the product name.

- ② Sample weight entry display

Weight



Piece Weight



Flashing weight value of the sample weight

in the unit (only g or oz are possible as units here)

-

_

Count
MEM 12
Adopt the value by briefly pressing the $On/Zero Off$ (Yes) key. Continue with point ③. Press the (Back) key to return to point ①.
③ Lower checking limit entry display Weight
Piece Weight
Check Weighing has been initialized. Piece Weight
Flashing quantity, if Check Counting has been initialized.
Pressing and holding the $Print Unit$ key will switch the unit between a weight unit and a pieces unit, if the product had not been initialized in advance (= empty product). This selection also affects point ④ (specifying upper limit).
Adopt the entered value by briefly pressing the On/Zero (Yes) key. Continue with point ④.
Press the (Back) key to return to point ②.
④ Upper checking limit entry display Weight

5E T.H I	
----------	--

Piece Weight

ig
angle flashing weight value of the upper limit, if

Check Weighing has been initialized.

Piece Weight

Flashing quantity, if Check Counting has

been initialized



For handling of an empty product, see point ③.

Adopt the value by briefly pressing the $O_{Off}^{On/Zero}$ (Yes) key. Continue with point (5).

Press the Target

 \int (Back) key to return to point ③.

The limit values are checked under the same conditions as defined in Chapter 5.2.4.2.1.

If a condition is not complied with, [LIM.ERR] is displayed briefly on the "PW" display.



Then the entry process is continued for the lower limit with the old value.

- ⑤ Saving entry display

Weight

<u>\| [</u> M. flashes

The "PW" display is empty.

Count



Adopt the product by briefly pressing the	On/Zero Off	(Yes) key or	Mem
raope are produce by briefly pressing are		(100) hey of	

Press the (Back) key to return to point (4).

- The applications are initialized depending on the product data. If Check Weighing has
 previously been deactivated and the product now contains a limit value > 0, Check
 Weighing is automatically initialized and vice versa.
- The applications are run, the change or initialization process has been completed.

5.2.8.4 Deleting product memory/parameters

- 1. Enter a product number via the keypad or using a scanner.
 - ▷ The number is displayed in the "Count" display.
- 2. Briefly press the key.
 - ▷ The following is displayed in the displays:



- 3. Briefly press the $\begin{bmatrix} C \\ +/- \end{bmatrix}$ key.
 - ▷ The following is displayed in the displays:

Weight



The "PW" display is empty.

4. Briefly pressing the On/Zero Off (Yes) key or the Mem (Exit) key will delete the product memory; briefly pressing the M+Menu (Exit) key will cancel the process.

5.2.8.5 Printouts

In the menu: [PRINT]- [CONTNT] it is possible to specify that the product memory and product name need to be printed:


Printout		Description
PROD-ID:	2	Product memory number (ID), here 2
PROD-NAME:	Metal4712123	Product name; max. 12 alphanumeric characters.

5.2.9 Handling of the second scale

4

Σ

If a second scale is activated, the "Weight" display looks like this: Weight



 $\overline{\Delta \lambda} 2$ scale WP2 is activated

Press and hold the $\begin{bmatrix} * \\ 1 & 2 \end{bmatrix}$ key to switch between the scales: If scale 1 is activated, this switches to scale 2 (and vice versa).

The scale values of the active scale are adopted into the display and used for the applications Counting, Checking as well as Printing and PC Output. This means that only the weights from the visible scale are used.

The initialization of the Counting application only takes place on scale 1. The initialization of the checking application, on the other hand, can take place on both scales. The initialized application is run on scale 2. A product can be loaded irrespective of the active scale.

5.3 Menu

Activated menu selections are marked with the following icon: [o].

When the last menu level has been reached, the active selected parameter, which is marked with [o], is displayed first.

If there has been a change in one or more menu items, these parameters are retained after exiting the menu.

If the menu is activated by pressing and holding the M^+_{Menu} key, then [M.E.N.U] is displayed. If the key is released, the first element of the top menu level [METRO] is displayed.

The menu item [END] can be used to exit a menu level. The display switches to the level above. If the top menu level is active and [E.N.D] is selected, the menu is exited.

Entries via the keyboard are only possible in a status in which an entry is expected.

The scale settings can be adjusted in the user menu (menu mode).

Note:

If appropriate interface options are installed, additional sub-menus may be available. Information on this can be found in the manual for the interface used.

5.3.1 Menu display in displays


Display of the first three menu levels in the three displays.

Weight	e menu levels in	Piece Weight
First menu lev	/el	Second menu level
Count)
Third menu	level	
Example:		J
- First menu level:		
Weight		First menu level
"PW" display and "(Count" display a	
 The second menu let 		e empty.
Weight		
MET		
Piece Weight		First menu level
].P.A	
		Second menu level
The "Count" display		
Weight	acheu.	
	рп	
	IN [_]	First an and beaut
Piece Weight		First menu level
	3.2.2	
		Second menu level
Count		
	0	

_____ Third menu level = selection

If there are more than three menu levels, the last level is displayed.

Second menu level	Third menu leve
Count	
Fourth menu level	
Example:	
- The second level is reached:	
Weight	
PP	
	Second menu level
Piece Weight	
	Third menu level
The "Count" display is empty.	
- The third level is reached:	
Weight	
	Second menu level
Piece Weight	
	Third menu level
The "Count" display is empty.	
- The fourth level is reached:	
Weight	
	Second menu level
Piece Weight	
	Third menu level



5.3.2 Menu mode

Call up the menu mode:

- 1. Press and hold the M+ Menu key.
 - In the "Weight" display, [M.E.N.U] is displayed for two seconds and then the 1st menu item [Metro] is displayed.

Weight

	M.E.N.LI
Weight	METRO

- 2. Press the Print Unit (No) key to access the next menu item.
 - ▷ The second menu item [UNIT] is shown in the "Weight" display.

Weight

veignt				
	 	T I		

- 3. Or press the Or/Zero (Yes) key to display the sub-menu in the "PW" display (in this example [[METRO][STAB.RA]]).
 - ▷ The second menu level is displayed in the "PW" display.

Weight	Piece Weight
METRO	

The "Count" display is empty.

- 4. To call up a menu selection (in this example [[METRO][STAB.RA]]), press the On/Zero Off (Yes) key.
 - ▷ The value 1D flashes in the "Count" display.

	Weight Piece Weight
	Count
	The current selection is marked with [o].
5.	Press the Unit (No) key in order to change the setting, or the Off (Yes) key in order to adopt the setting.
	The next element of the second menu level [FILTER] is displayed in the "PW" display.
	Weight Piece Weight
	The "Count" display is empty.
6.	Press the Print Unit (No) key in order to select the next sub-menu in the "PW" display, or press the (Back) key to go back.
	Weight Piece Weight
7.	Press the On/Zero Off Off (Yes) key when [END] is displayed. ▷ The second menu item [UNIT] is shown in the "Weight" display.
	Weight
8.	Press the Menu (Exit) key to exit the menu.
	tered values in the menu
-	For more than the second s

For menu items with numerical settings, e.g. interval, the current setting is displayed with flashing digits.

1. Press the On/Zero Off (Yes) key in order to adopt the setting, or the Unit (No) key in order to continue processing.

5.3.3

- 2. To exit the menu elements, press the Or/Zero Off (Yes) key to access the next menu item, **Print** Unit
 - (No) key to access the top level of the current menu. or the

5.3.4 **Menu navigation**

Overview of the first menu level:

– METRO	Metrology (see Chapter 5.3.4.1)
– UNIT	Weight units (see Chapter <mark>5.3.4.2</mark>)
– OP.FUNC	Operating functions (see Chapter 5.3.4.3)
– PRINT	Printer outputs (see Chapter 5.3.4.4)
– PRN.COM	Printer port communication (see
	Chapter
– PC.OUT	PC output (see Chapter <mark>5.3.4.6</mark>)
– PC.COM	PC port communication (see Chapter 5.3.4.7)
— CAL.ADJ	Calibration/adjustment (see Chapter 5.3.4.8)
– AD.CON2	ADC configuration of the second scale, is only
	displayed if [AD.CON2] (see Chapter 5.3.4.9)
	has been activated previously
— INFO	Info (display of serial number and type
	designation)(see Chapter <mark>5.3.4.10</mark>)
– SECURE	Block menu items (see Chapter 5.3.4.11)
– E.N.D.	Exiting menus

[METRO] menu selection 5.3.4.1

The functions of the displays and scales can be adjusted in this menu. Factory settings are marked with "*"

METRO	
– STAB.RA	Stability range value for both scales
– 0.5D	1/2 d
– 1D	1 d*
— 2D	2 d
– 4D	4 d
– FILTER	Adjustment filter valid for both scales
LOW	Lower accuracy, short stabilization time
— MED	Normal accuracy, average stabilization time*
– ні	High accuracy, long stabilization time
– A.ŻERO.T	Automatic zero point tracking valid for both
	scales
– OFF	Switching off
— 0.5D	Drift up to 0.5 d/second*
– 1D	Drift up to 1 d/second
— 3D	Drift up to 3 d/second
– AUT.OFF	Counter for automatic switching off
– OFF	Switching off*
— 1 MIN	Switching off after 1 minute with no activity
— 5 MIN	Switching off after 5 minutes with no activity
— 10 MIN	Switching off after 10 minutes with no activity
– RESET	Factory settings
— NO	not reset*
– YES	activated
- END	Exit menu level

5.3.4.2 [UNIT] menu selection

The weight unit can be selected in this menu. Factory settings are marked with "*"

UNIT	
— kg	Kilogram
- OFF	deactivated
– ON	activated*
— g	Gram
- OFF	deactivated
– ON	activated*
– Ib '	Pound
– OFF	deactivated
– ON	activated*
— oz	Ounce
- OFF	deactivated
– ON	activated*
— lb:oz	Pound:ounce
- OFF	deactivated*
– ON	activated
– RESET	Factory settings
- NO	not reset*
– YES	activated
– END	Exit menu level

5.3.4.3 [OP.FUNC] menu selection

The scale parameters can be specified in this menu.

Factory settings are marked with "*"

OP.FUNC		
— WP2	Activate a second scale	
– OFF	deactivated*	
– ON	activate	
– UNIT.ON	Select unit when switching on	
	Last unit used when switching off*	
— kg	Kilogram	
— g	Gram	
— Ib	Pound	
— oz	Ounce	
— lb:oz	Pound:ounce	
– ZERO.RA	Zero range	
– 2%	2% max. load	
— 10%	10% max. load*	
– A.OPT	Automatic optimization of the sample weight	
– OFF	deactivated	
— ON	activate*	
– A.TARE	Automatic tare	
– OFF	deactivated*	
– ON	1st stable weight is tared	
– ON-ACC	Stable loads within the tolerance limits are	
	tared (in all checking applications)	
– BEEP.OP		

⊢ OFF ⊢ ON − BEEP.SI ⊢ OFF − ACCEPT	Signal when the sample weight has been automatically optimized deactivated activate* Signal (in the Checking application) deactivated* Signal when the weight is within the tolerance limits
– UNDER	Signal when the weight is below the lower limit
– OVER	Signal when the weight is above the upper limit
– UNDOVR	Signal when the weight is outside the tolerance
	limits
- BEEP.KE	Key tone
	deactivated
	activated*
	Totalizing setting deactivated
– OFF – AUTO	
	Automatic totalizing Manual totalizing*
	Duration of the background lighting (D.LIGHT =
	AUTO)
⊢ 3 SEC	Switching off of the background lighting after 3
	seconds with no activity
— 5 SEC	Switching off of the background lighting after 5
	seconds with no activity*
— 8 SEC	Switching off of the background lighting after 8
	seconds with no activity
– D.LIGHT	Background lighting of the display
- OFF	deactivated
	activated
– AUTO	Switches on when a key is pressed or the
	displayed weight changes* Communication module
⊢ COM.EQU └─ OFF	deactivated*
BLUE.TH	Bluetooth activated (if the Bluetooth module is
	installed)
	WiFi activated (if the WiFi module is installed)
– ETHER.N	Ethernet activated (if the Ethernet module is
	installed)
– RESET	Factory settings
- NO	not reset*
– YES	activated
- END	Exit menu level

5.3.4.4 [PRINT] menu selection

print conditions and printouts can be configured in this menu.

Factory settings are marked with "*"

PRINT	
– STABLE	Print criteria
⊢ OFF	Values are printed immediately
– ON	Values are only printed if they are stable*

- A.PRINT	Automatic printout
– OFF	deactivated*
– ON.STAB	One-time printing on stability, if weight > 1 d
— INTER	Printing in the specified interval
- 13600	1 3,600 seconds
– CONT	Print continuously
– ACCEPT	Printing on stability and within the tolerance
	limits
- CONTNT	Content of a printout
– RESULT	Displayed value
– OFF	deactivated
— ON	activated*
– GROSS	Gross value
– OFF	deactivated*
	activated
- NET	Net value
- OFF	deactivated*
	activated
	Tare
- OFF	deactivated*
	activated
- PRE.TAR	Preset tare value, if available
- OFF	deactivated*
	activated
— HEADER	Header is described via the PC interface (see
	SBI specifications) deactivated*
	activated
│	
FOOTER	Info on the footer; is specified via the user interface (see SBI specifications)
	deactivated*
	activated
– MODE	Info on the application mode
	deactivated*
- ON	activated
	Application parameters
⊢ OFF	deactivated*
– ON	activated
– PW	Sample weight
- OFF	deactivated*
— ON	activated
– TOTAL	Totalizing memory/statistics data
- OFF	deactivated*
– RESULT	Totalizing memory is printed
	Result and statistics parameters such as min,
	max, etc.
– PROD.ID	Product memory location number
- OFF	deactivated*
	activated
– PROD.NA	Product name; max. 12 alphanumeric
	characters
- OFF	deactivated*

- SERNO - SERNO - OFF - ON - LI.SET - FORMAT - MULTI - SINGLE - FEED - LINE - 4LF - FORM - END	activated Serial number of the active (displayed) scale Example: W1 38457989 or W2 38457989 deactivated* activated Totalizing memory Format sent to printer and PC Multi-line (single-column) printout* Single-line printout; the entire content defined above is printed in one line Setting of the paper feed One-line feed Four-line feed* Page feed after printing
– PROD.LI	Print entire product memory
	deactivated*
— YES	Print all products
	This happens at the time when Yes is selected. Then all products are printed. Afterwards No is specified automatically.
– RESET	Factory settings
- NO	not reset*
- YES	activated
⊢ END	Exit menu level

In this case "without stability" means weight values with or without stability. In this case "with stability" means only weight values with stability.

5.3.4.5 [PRN.COM] menu selection

The parameters for the print communication can be specified in this menu. Factory settings are marked with "*"

PRN.COM	
– BAUD	Baud rate
- 2400	2,400
- 4800	4,800
- 9600	9,600*
- 19200	19,200
- 38400	38,400
- 57600	57,600
- 115200	115,200
– PARITY	Parity
— 7 EVEN	7 data bits, even parity
— 7 Odd	7 data bits, odd parity
– 7 NONE	7 data bits, no parity
– 8 NONE	8 data bits, no parity*
– STOP	Stop bit
	1*
- 2	2
– RESET	Factory settings
- NO	not reset*
– YES	activated

5.3.4.6

	- END	Exit menu level	
[PC.OUT] menu selection			
	The parameters for the PC output can be	specified in this menu.	
	This is an additional interface based on R	S-232 with USB-C port.	
	Factory settings are marked with "*"	·	
	PC.OUT		
	– MODE	PC output mode	
	- OFF	deactivated*	
	– MAN.OUT	Manual output without stability = device	
		expects a request via the serial interface.	
		The response is generated immediately .	
	— MAN.STA	Manual output without stability = device	
		expects a request via the serial interface.	
	– INT.OUT	Interval output with 16 characters without	
		stability = device sends weight values	
		automatically without stability in the selected	
		interval of the update cycles for the display (see	
		below)	
	– AUT.OUT	Automatic output without stability with 16	
		characters = device sends weight values	
		without stability automatically on every update	
		cycle for the display	
	– AUT.STA	Automatic output with stability with 16	
		characters on every value change = device	
		automatically sends weight values with stability	
		Specify output interval (if INT.OUT is selected)	
	– 1 CYC – 2 CYC	Every display cycle After 2 display cycles	
		After 5 display cycles	
		After 10 display cycles*	
	- 20 CYC	After 20 display cycles	
	- 50 CYC	After 50 display cycles	
	– 100 CYC	After 100 display cycles	
	– RESET	Factory settings	
	⊢ NO	not reset*	
	– YES	activated	
	— END	Exit menu level	

In this case "without stability" means weight values with or without stability. In this case "with stability" means only weight values with stability!

5.3.4.7 [PC.COM] menu selection

The parameters for the PC communication can be specified in this menu. Factory settings are marked with "*"

Baud rate
4,800
9,600*
19,200
38,400

- 57600	57,600
- 115200	115,200
– PARITY	Parity
— 7 EVEN	7 data bits, even parity
— 7 Odd	7 data bits, odd parity
— 7 NONE	7 data bits, no parity
– 8 NONE	8 data bits, no parity*
– 7 MARK	7 data bits, mark parity
– 7 SPACE	7 data bits, space parity
– STOP	Stop bit
	1*
- 2	2
– HAND.SH	Handshake
- NONE	No handshake*
— XON.XOF	No function
– RESET	Factory settings
— NO	not reset*
– YES	activated
— END	Exit menu level
-	

5.3.4.8 [CAL.ADJ] menu selection

Factory settings are marked with "*"

CAL.ADJ	
– CAL	Initiates a two-point calibration (zero and max. load)
- WP 1	Initiates a two-point calibration for WP1
	The calibration weight can be freely selected by entering the value.
— WP 2	Only visible if a second scale is defined in the menu [OP.FUNC]- [WP2]:
	Initiates a two-point calibration for WP2
	The calibration weight can be freely selected by
	entering the value.
	Initiates a linearization (zero, half weighing range and max. load)
- WP 1	Initiates a linearization for WP1
	The weight values are fixed.
— WP 2	Only visible if a second scale is defined in the
	menu [OP.FUNC]- [WP2]:
	Initiates a linearization for WP2
	The weight values are fixed.
— GEO	The adjustment of the calibration based on the current location is carried out using the
	geographic adjustment factor (GEO)
- 12	Selectable range from 0 to 31 in individual increments*
	Exit menu level

5.3.4.9 [AD.CON2] menu selection

The ADC configuration for scale 2 is carried out in this menu.

Note:

- The increment (readability) of the second scale must be bigger than the increment of the first scale (d2 > d1).

Example: If the readability of the first scale = d1 = 1 g, then the readability of the second scale must be = d2 = 2 g (d2 cannot be 1 g or 0.5 g).

- The quotient from calibration weight and readability (CAL.WGT / d) must be \geq 1,000. This also affect the max. load (MAX.CAP), as: MAX.CAP \geq CAL.WGT.

AD.CON2	
⊢ D	The readability (1 position) can only be changed
	Print Target
	via the \bigcup_{nit} key or the \bigcup_{nit} key.
- CAL.WGT	Calibration weight, which is used in [CAL.ADJ]-
	[CAL]- [WP2] and can be freely selected within
	30% of the max. load and min. load
- MAX.CAP	Max. load of the scale
	If this entered load has been reached, "H" is
	displayed.
- END	Menüebene verlassen
5.3.4.10 [INFO] menu selection	
INFO	
– SER.NUM	Display serial number (if two scales are
	selected, two serial numbers are displayed)
⊢ W1.471	Print Target
	Using the Unit (No) key or the key
	will display the next six digits or display the first
	six digits again.
— W2.471	Using the Unit (No) key or the key
	will display the next six digits or display the first
— END	six digits again.
	Display type name of the scale
- SFT.VER	Display type name of the scale
- SF1.VEK	Display software version
— GEO.DAT	Display selected GEO area, which is valid for
	both scales
⊢ 12	Selected GEO area
	Info WP1
	Readability
	Max. load
– END	
— WP2	Info WP2, if selected in the menu
	Accuracy
	Max. load
- END	Max Iouu
– LFT.MOD	Display status of the legal metrology mode: ON
	or OFF for both scales
	Default mode
	- craat mode

⊢ ON – LOCK.SW	Legal metrology mode is activated Display status of the locking switch for both
	scales
- OPEN	Switch is open
– OPEN – CLOSED	Switch is closed
– END	Exit menu level

5.3.4.11 [SECURE] menu selection

The safety setting (lock) for menu access can be defined via this menu in order to prevent unauthorized interventions.

Factory settings are marked with "*"

SECURE	
– S.METRO	Metrology menu
- OFF	Menu item can be changed*
– ON	Lock menu item
– S.UNIT	Unit menu
- OFF	Menu item can be changed*
– ON	Lock menu item
– S.OP.FUN	Operating functions menu
- OFF	Menu item can be changed*
– ON	Lock menu item
— S.PRINT	Print menu
– OFF	Menu item can be changed*
— ON	Lock menu item
– S.PR.COM	Print communication menu
– OFF	Menu item can be changed*
— ON	Lock menu item
– S.PC.OUT	PC output menu
– OFF	Menu item can be changed*
– ON	Lock menu item
– S.PC.COM	PC communication menu
– OFF	Menu item can be changed*
- ON	Lock menu item
– S.CAL.AD	Calibration menu
– OFF	Menu item can be changed*
- ON	Lock menu item
– S.ADC.CO	ADC menu
– OFF	Menu item can be changed*
- ON	Lock menu item
- RESET	Restore factory setting of the current menu
– OFF	Menu item can be changed*
	Lock menu item
- END	Exit menu level


5.4 Calibration, adjustment

5.4.1 Calibration

The calibration takes place following activation of menu item [WP1] or [WP2] in the menu. Weights for calibration points can be freely selected when the value is flashing. If the unit used for normal weighing is pounds, the calibration unit is also lb (pound). If the unit used for normal weighing is ounces, the calibration unit is also oz (ounce). If a metric unit (kg or g) is used for normal weighing, the calibration unit is kg (kilograms).

5.4.1.1 Calibration

- 1. WP1 or WP2 (if available) has been activated.
 - ▷ The following is displayed in the displays:



- 2. Remove all weights from the weighing pan and press the off (Yes) key.
 - \triangleright Now the first calibration point is adopted and saved internally.


In this time, the following is displayed in the displays:

Weight




Second menu level = calibration

Piece Weight



Third menu level = here WP1

Count



ig J Wait message, until the weight has been


adopted

The second weighing point is displayed flashing.





Second menu level = calibration



via the keypad

3. If the correct second weighing point is selected, press the On/Zero Off (Yes) key.

Adopt the second weighing point. During the internal calculations, the same display appears in the display as for step 2.

If the calculations were error-free, the following is displayed on the display for around 2 seconds:

Weight

 Second menu level = calibration

Piece Weight



Third menu level = here WP1

Count



OK message = calibration is complete

The device exits the menu and returns back to the normal application mode.

If an error occurred during the calibration, the following is displayed on the displays for around 2 seconds:

Weight	
	Second menu level = calibration
Piece Weight	
	Third menu level = here WP1
Count	
[ALERR	

Error message = calibration is canceled

The device exits the menu and returns back to the normal application mode.

The process is repeated for WP2. WP2 is displayed in the "PW" display and the device switches automatically to scale 2.

5.4.1.2 Menu settings

CAL.ADJ	
– CAL	Selection of the scale to be calibrated and start of calibration
— WP 1	Initiates a two-point calibration (zero and max. load) for WP1
	This process can be canceled by pressing the "Exit" key.
	The calibration weight can be freely selected by entering the value.
— WP 2	Only visible if a second scale is defined in the
	menu:
	Initiates a two-point calibration (zero and max. load) for WP2
	This process can be canceled by pressing the "Exit" key.
	The calibration weight can be freely selected by entering the value.
– LIŃ	-
– GEO	
- END	Exit menu level

5.4.2 Linearization

The linearization takes place following activation of menu item [WP1] or [WP2] in the menu.


5.4.2.1 Linearization process

Linearization weights are fixed. There are three points: Zero, 0.5* maximum load, maximum load.

The procedure is similar to the procedure for the calibration. The following is displayed in the "Weight" display: [LIN]


- ► WP1 or WP2 (if available) has been activated.
 - ▷ The following is displayed in the displays:

Weight



Third menu level = here WP1

5.4.2.2



5.4.3 Geographic adjustment factor (GEO) – procedure

The adjustment of the calibration based on the current location is carried out using the geographic adjustment factor [GEO]. (Settings from 0 ... 31 are available.) The table under 5.4.4 contains the GEO values for a wide range of latitudes.

Call up the menu mode:

- 1. Press and hold the Menu key.
 - ▷ The first menu item [METRO] is displayed in the "Weight" display.
 - Weight



Press the Print Unit (No) key until the menu item [CAL.ADJ] is displayed in the "Weight" display.

Weight



3. Press the On/Zero (Yes) key to access the sub-menu item.

Print Unit (No) key until the menu item [GEO] is displayed in the "PW" display. 4. Press the



- On/Zero Off (Yes) key to start the GEO selection. 5. Press the
 - ▷ The GEO value [12] set by default flashes in the display.



6. To change the value, a value from 0 ... 31 can be entered using the keypad.

- 7. Press the On/Zero Off Off (Yes) key to confirm the GEO value.
 - ▷ The GEO value has been saved when [END] is displayed in the "PW" display.

Weight	Piece Weight

- On/Zero 8. Press the Off (Yes) key.
 - The menu item [INFO] is displayed in the "Weight" display. \triangleright
- 9. Press the (Exit) key to exit the menu.

GEO code table 5.4.4

						Altit	ude in m	eters				
		0	325	650	975	1,300	1,625	1,950	2,275	2,600	2,925	3,250
		325	650	975	1,300	1,625	1,950	2,275	2,600	2,925	3,250	3,575
						Alt	itude in f	feet				
		0	1,016	2,130	3,200	4,260	5,330	6,400	7,460	8,530	9,600	10,660
		1,060	2,130	3,200	4,260	5,330	6,400	7,460	8,530	9,600	10,660	11,730
Lat	itude					(GEO valu	e				
0°00'	5°46'	5	4	4	3	3	2	2	1	1	0	0
5°46'	9°52'	5	5	4	4	3	3	2	2	1	1	0
9°52'	12°44'	6	5	5	4	4	3	3	2	2	1	1

						Altit	ude in m	eters				
		0	325	650	975	1,300	1,625	1,950	2,275	2,600	2,925	3,250
		325	650	975	1,300	1,625	1,950	2,275	2,600	2,925	3,250	3,575
						Alti	itude in f	feet				
		0	1,016	2,130	3,200	4,260	5,330	6,400	7,460	8,530	9,600	10,660
		1,060	2,130	3,200	4,260	5,330	6,400	7,460	8,530	9,600	10,660	11,730
Lati	tude					(GEO valu	е				
12°44'	15°06'	6	6	5	5	4	4	3	3	2	2	1
15°06'	17°10'	7	6	6	5	5	4	4	3	3	2	2
17°10'	19°02'	7	7	6	6	5	5	4	4	3	3	2
19°02'	20°45'	8	7	7	6	6	5	5	4	4	3	3
20°45'	22°22'	8	8	7	7	6	6	5	5	4	4	3
22°22'	23°54'	9	8	8	7	7	6	6	5	5	4	4
23°54'	25°21'	9	9	8	8	7	7	6	6	5	5	4
25°21'	26°45'	10	9	9	8	8	7	7	6	6	5	5
26°45'	28°06'	10	10	9	9	8	8	7	7	6	6	5
28°06'	29°25'	11	10	10	9	9	8	8	7	7	6	6
29°25'	30°41'	11	11	10	10	9	9	8	8	7	7	6
30°41'	31°56'	12	11	11	10	10	9	9	8	8	7	7
31°56'	33°09'	12	12	11	11	10	10	9	9	8	8	7
33°09'	34°21'	13	12	12	11	11	10	10	9	9	8	8
34°21'	35°31'	13	13	12	12	11	11	10	10	9	9	8
35°31'	36°41'	14	13	13	12	12	11	11	10	10	9	9
36°41'	37°50'	14	14	13	13	12	12	11	11	10	10	9
37°50'	38°58'	15	14	14	13	13	12	12	11	11	10	10
38°58'	40°05'	15	15	14	14	13	13	12	12	11	11	10
40°05'	41°12'	16	15	15	14	14	13	13	12	12	11	11
41°12'	42°19'	16	16	15	15	14	14	13	13	12	12	11
42°19'	43°26'	17	16	16	15	15	14	14	13	13	12	12
43°26'	44°32'	17	17	16	16	15	15	14	14	13	13	12
44°32'	45°38'	18	17	17	16	16	15	15	14	14	13	13
45°38'	46°45'	18	18	17	17	16	16	15	15	14	14	13
46°45'	47°51'	19	18	18	17	17	16	16	15	15	14	14
47°51'	48°58'	19	19	18	18	17	17	16	16	15	15	14
48°58'	50°16'	20	19	19	18	18	17	17	16	16	15	15
50°16'	51°13'	20	20	19	19	18	18	17	17	16	16	15

						Altit	ude in m	eters				
		0	325	650	975	1,300	1,625	1,950	2,275	2,600	2,925	3,250
		325	650	975	1,300	1,625	1,950	2,275	2,600	2,925	3,250	3,575
						Alt	itude in f	feet				
		0	1,016	2,130	3,200	4,260	5,330	6,400	7,460	8,530	9,600	10,660
		1,060	2,130	3,200	4,260	5,330	6,400	7,460	8,530	9,600	10,660	11,730
Lati	tude					C	GEO valu	e				
51°13'	52°22'	21	20	20	19	19	18	18	17	17	16	16
52°22'	53°31'	21	21	20	20	19	19	18	18	17	17	16
53°31'	54°41'	22	21	21	20	20	19	19	18	18	17	17
54°41'	55°52'	22	22	21	21	20	20	19	19	18	18	17
55°52'	57°04'	23	22	22	21	21	20	20	19	19	18	18
57°04'	58°17'	23	23	22	22	21	21	20	20	19	19	18
58°17'	59°32'	24	23	23	22	22	21	21	20	20	19	19
58°17'	59°32'	24	23	23	22	22	21	21	20	20	19	19
60°49'	62°90'	25	24	24	23	23	22	22	21	21	20	20
62°90'	63°30'	25	25	24	24	23	23	22	22	21	21	20
63°30'	64°55'	26	25	25	24	24	23	23	22	22	21	21
6 4°55'	66°24'	26	26	25	25	24	24	23	23	22	22	21
66°24'	67°57'	27	26	26	25	25	24	24	23	23	22	22
67°57'	69°35'	27	27	26	26	25	25	24	24	23	23	22
69°35'	71°21'	28	27	27	26	26	25	25	24	24	23	23
71°21'	73°16'	28	28	27	27	26	26	25	25	24	24	23
73°16'	75°24'	29	28	28	27	27	26	26	25	25	24	24
75°24'	77°55'	29	29	28	28	27	27	26	26	25	25	24
77°55'	80°56'	30	29	29	28	28	27	27	26	26	25	25
80°56'	85°45'	30	30	29	29	28	28	27	27	26	26	25
85°45'	90°00'	31	30	30	29	29	28	28	27	27	26	26

5.5 Print

During every print process, [PRINT] is displayed in the "Weight" display for around 1 second (maximum).

Weight

$$PR + NT$$

5.5.1 Output format

General definitions:

- Space = 0x20 hex as ASCII character.
- CRLF = 0x0D0A hex as ASCII character.
- The algebraic sign is inserted directly next to the MSB of a value with a space in between.

```
Example: - 14.112
```

- The values are written over each other:

Example:

```
    - 10.075 oz OVERCRLF
    0.015 oz TCRLF
    - 10.060 oz G#CRLF
```

- A stable weight value is printed with a unit. 10.25 kg

A non-stable weight value is printed without a unit: 10.25

- A stable negative gross weight value has a "!" instead of a unit: - 0.25 !

A non-stable negative gross weight value has no unit symbol: - 0.25

- A defective format is generated if no valid weight value is displayed/if an error message is displayed. The error number in the display is also printed: e.g. [Err 8.1].

5.5.2 Product memory

Print entire product memory, see menu item [PRINT]- [PROD.LI]:

Printout	Description
PROD-ID: 2	Product memory, here 2. Memory 1 is empty
PROD-NAME: Metal4712123	Product description
5.23 g PT	Preset tare value
WREF 0.43245 g	Sample weight
UNDER LIMIT 580.0 g	Lower limit
OVER LIMIT 681.0 g	Upper limit
MODE: CHECK WEIGHT	Checking application
	Line feed
PROD-ID: 3	Product memory, here 3. Only counting without preset tare value
PROD-NAME: Metal458ab	Product description
0.0 g PT	Preset tare value
WREF 0.83241 g	Sample weight
UNDER LIMIT 0.0 g	Lower limit
OVER LIMIT 0.0 g	Upper limit
MODE: CHECK OFF	Checking application
	Line feed
PROD-ID: 8	Product memory, here 8. Memory 4, 5, 6, 7 are empty
PROD-NAME: Metal471	Product description
0.0 g PT	Preset tare value
WREF 0.0 g	Sample weight
UNDER LIMIT 580.0 g	Lower limit
OVER LIMIT 681.0 g	Upper limit
MODE: CHECK WEIGHT	Checking application
	Line feed

5.5.3 Menu settings: Print

The [PRINT] menu item has the following sub-menus:

PRINT	
– STABLE	Print criteria
– A.PRINT	Automatic printout
– CONTNT	Content of a log printout
– LI.SET	Frame format
– PROD.LI	Print entire product memory
– RESET	Set [PRINT] menu item to factory setting
– END	Exit menu level

Print conditions

Define conditions in the [PRINT] menu item (see Chapter 5.3.4.4)

Printout contents

- The contents of a print log are specified in the menu item [PRINT]- [CONTNT] (see Chapter 5.3.4.4).

Any parameter selected here is printed for **every** printout, even if there is no content there.

Example: If there is no tare value available, 0.000 kg T is printed.

This does not apply for totalizing and a preset tare value. The parameters for totalizing are only printed if totalizing is activated in the menu!

The preset tare is only printed if the memory is published.

 In the case of a stable weight value, the unit is always printed. A selection is not possible!

Frame format

The frame format is configured in the menu item [PRINT]- [LI.SET] (see Chapter 5.3.4.4).

Settings for the print interface

Selection of the interface parameters for printing [PRN.COM] (see Chapter 5.3.4.5).

5.6 PC output

This is an additional interface, which is based on the USB-C standard.

5.6.1 SBI interface

A computer connected via the PC interface (SBI communication) can send control commands to the analysis device in order to control the scale or application functions. All commands have a shared frame format (data input format). They start with the characters ESC and end with the command end EOC (end of command). The end of command may also be a combination of CR and LF. The scale ignores all entries after EOC and before ESC.

Reading the displayed value:

ESC	-				Р					EOC					-
Res	ponse	(16 by	/tes):												_
٧	W	W	W	W	W	W	W	W	W	E	E	E	CR	LF	_

V	Algebraic sign	Possible characters: "+", "–", " "
W	Weight value	Possible characters: "0""9", ". ", " "
E	Unit	Possible characters: "a""z", "A""Z", " "
CR	Carriage return	ASCII 0x0D
LF	Line feed	ASCII 0x0A

This format is also used for automatically generated telegrams, which are released according to the menu settings: [INT.OUT], [AUT.OUT], [AUT.STA] (see above).

Zeroing the scale:

ESC Z	EOC
-------	-----

Response: see special response telegrams

Taring the scale:

	ESC	Т	EOC	
--	-----	---	-----	--

Response: see special response telegrams

Special response telegrams:

There are some special responses, which are used as standard responses. Example: Error or confirmation. Special response telegrams are always 5 bytes.

OK (confirmed)

1	2	3	4	5
0	К	!	CR	LF

The scale confirms error-free performance of the command.

ERROR (error)

1	2	3	4	5
E	R	R	CR	LF

The scale reports an error when performing the command.

LOCKED (locked)

1	2	3	4	5
L	0	С	CR	LF

The command cannot be performed because a parameter is currently blocked.

6 Maintenance/repairs/cleaning

6.1 Repairs

Disconnect a defective device from the mains immediately.

Defective or damaged cables or screw connections must be replaced as a complete unit.

▲ WARNING

Improper repairs can pose considerable risks to the user.

Only have repairs carried out by Minebea Intec qualified dealers using original spare parts.

6.2 Cleaning

6.2.1 Instructions for cleaning

The device must be cleaned of contaminants on a regular basis.

Before cleaning, maintenance, or repairs, disconnect the device from the supply voltage.

If the scale is in a dry environment, wipe the weighing platform with a damp cloth. Household cleaning agents can be used. Please check the information provided by the manufacturer.

In the case of devices with an IP43 protection grade, no liquid must get into the scale.

The device must not be cleaned using a high-pressure or steam cleaner. Observe the IP protection grade.

If the device is cleaned with water that is too hot or too cold due to temperature differences, condensation may form in the device. Condensation may cause malfunctions in the device.

6.2.2 Cleaning agents

NOTICE

Some cleaning agents may not be compatible with the device material.

- Only use disinfectants and cleaning agents in line with the manufacturer's instructions.
- Do not use cleaning agents that are very acidic, very alkaline, or that contain a high level of chlorine. Avoid substances with a high or low pH value as otherwise there is an increased risk of corrosion.
- Do not use any abrasive sponges containing iron, steel brushes, or cleaning sponges made of steel wool.
- Always test cleaning agents and materials in non-critical areas first before using them.

7 Waste disposal policy

If the packaging is no longer required, please take it to your local waste disposal facility and/or a reputable disposal company or collection point. The packaging largely consists of environmentally friendly materials, which are suitable for recycling.

It is not permitted—even for small businesses—to dispose of this product with the regular household waste or at collection points run by local public waste disposal companies.

EU legislation requires its Member States to collect electrical and electronic equipment and dispose of it separately from other unsorted municipal waste so that it can then be recycled.

Before disposing of or scrapping the product, any batteries should be removed and taken to a suitable collection point.

Please see our T&Cs for further information.

We reserve the right not to accept products that have been contaminated with hazardous substances (ABC contamination) for repair.

8 Error correction

The table lists frequent problems, as well as possible causes and corrective measures. If the problem persists, inform Minebea Intec or an authorized dealer.

Symptom	Possible cause	Corrective measure
Switching on not possible	Scale is not supplied wi- th power	Check connections and voltage
Poor accuracy	Incorrect calibration Unstable environment	Perform a calibration Put scale in a suitable location
Application cannot be called up	Application has not be- en activated	Activate the application in the menu
Unit cannot be cal- led up	Unit has not been acti- vated	Activate the unit in the menu
Battery icon is flas- hing	Low battery level	Connect scale to the mains and charge battery
[Err 8.1]	Error during switching on	Read weight exceeds start-up/zeroing limit
[Err 8.2]	Error during switching on	Read weight falls below start-up/ze- roing limit
[Err 8.3]	Overload range error	Read weight exceeds overload limit
[Err 8.4]	Underload range error	Read weight falls below overload limit
[Err 8.5]	Tare outside of the tare range	Adjust tare value accordingly
[Err 8.6]	Display capacity fallen short of	Weight > 6 characters
[Err 9.5]	Incorrect calibration da- ta	Repeat calibration
	Busy	Display during tare setting, zero point setting, printing
[NO]	Action not permitted	Function cannot be performed
[CAL.ERR]	Calibration error Unstable environment Incorrect calibration weight	Repeat calibration Put scale in a suitable location Use correct calibration weight
[LIM.ERR]	Entered value outside of the permitted range	Enter a permitted value
[REF.ERR]	Invalid reference weight	The weight on the load plate is too high or too low in order to define a valid re- ference weight. Reduce or increase re- ference weight
Battery cannot be fully charged	Battery is defective	Have battery replaced by authorized Minebea Intec service dealer.

8.1 Service information

Contact the authorized service partners if a problem cannot be rectified with the aid of the troubleshooting information or is not described there. Our website http://www.puroscales.com provides information about your closest service partner.

9 Technical data

9.1 Specification


Model number	EF - LT	P1	P3	P6	P15	P30
Max load (g)		1,500	3,000	6,000	15,000	30,000
Readability d (g)	- 30 d	0.05	0.1	0.2	0.5	1
Max. resolution	_	30,000	30,000	30,000	30,000	30,000
Readability d (g)	- 6 d	0.2	0.5	1	2	5
Max. resolution	-	7,500	6,000	6,000	7,500	6,000
Application package Counting Professional		Weighing, Automatic Tare, Automatic Printing, Counting Professio Checking				g Professional,
Min. recommended samp- le weight (g)		1	2	4	10	20
Min. average recommen- ded sample weight (g)	_	0.005	0.01	0.02	0.05	0.1
Max. internal resolution during counting		1:1,500,000	0 internal/1:30	0,000 external		
Weight units		kg, g, lb, oz	z, lb:oz			
Version/materials		Housing m 304	ade of ABS pla	stic, weighing p	olatform made	of stainless steel
Protection grade		IP43				
Display		3-window LCD display with white background lighting, digit height 1.1 inches/28 mm				digit height 1.1
Indicator displays		3 LEDs (yel signal	low, green, red), function can	be configured,	acoustic alarm
Memory		Library for	30 products			
Keypad		8 function,	12 mechanical	number keys		
Zero range		2 or 10% of	the max. load	of the scale		
Tare range		Max. load v	via subtraction			
Stabilization time		1 second				
Automatic zero point cor- rection		Off, 0.5; 1 c	or 3 display inci	rements		
Safe overload range		150% of the	e max. load of t	he scale		
Leveling aids		Externally	visible level ind	licator and adju	ıstable, non-sli	p leveling feet
Electrical supply		U _{DC} = 5 V, le lithium b		/60 Hz power	supply or insta	lled rechargeab-
Battery operation time		•	•	n time (with star urs charge time		between the
Calibration		External, w	ith freely selec	table calibratio	n weights	


Model number	EF - LT	P1	P3	P6	P15	P30		
Interface		USB-C, prir	USB-C, printer port, RS-232 installed, Bluetooth, or WIFI optional					
Operating temperature (°C)		-10 +40						
Storage temperature (°C)		-10 +50						
Product dimensions (W x D x H)		310 x 302 x	115 mm					
Platform size (W × L)	_	280 x 180 r	nm					
Shipping dimensions (W x D x H)		370 x 370 x	220 mm					
Net weight	_	2.9 kg						
Shipping weight	_	4.5 kg						

9.2 Accessories


Option	Order no.			
Data printer	YP-DP1			
Paper for data printer	YP-P1			
Printer cable	YP-CAS1			
USB-C cable	YP-CAC1			
Manual laser barcode scanner	MD2000			
Scanner cable	YP-CYSSR1			
USB charging device	YP-PS1			
Weighing hooks	YP-H1			

9.3 Dimensions






All dimensions in mm / inch



All dimensions in mm / inch



All dimensions in mm / inch

10 Appendix

10.1 Printouts

The elements to be printed are configured in the menu.

Weighing printout

Printout		Description
5.003 g	N	Positive net weight value
- 0.003 g	N	Negative net weight value
2.003 g	G	Positive measured gross weight value
2.003 g	G#	Positive calculated gross weight value
1.003 g	Т	Tare weight value (measured value)
0.010 g	PΤ	Preset tare value (entered value)
- 0.010 !	G	Stable gross weight below zero

Counting application printout (without checking application)

Printout	Description
441 pcs QNT	Positive value
- 41 pcs QNT	Negative value
MODE: COUNT	Activated application
WREF 4.15431 oz	Sample weight

Check Weighing application printout (Counting application not initialized)

Printout		Description
115	g OVER	Positive net value > upper limit
- 115	g OVER	Negative net value < negative upper limit
99	g ACCEPT	Positive net value in the target range
75	g UNDER	Positive net value < lower limit
MODE: CHECKW	EIGH	Activated application
UNDER LIMIT	81 g	Lower limit
OVER LIMIT	100 g	Upper limit

Check Weighing application printout (Counting application initialized)

Printout	Description
115 g OVE	Positive net value > upper limit
– 115 g OVE	R Negative net value < negative upper limit
99 g ACC	PEPT Positive net value in the target range
75 g UND	Positive net value < lower limit
115 pcs	Net number >
MODE: CHECKWEIGH	Activated application
UNDER LIMIT 81 g	Lower limit
OVER LIMIT 100 g	upper limit
WREF 0.35423 oz	Sample weight, as displayed

Check Weighing application printout with parts (Counting not initialized)

Printout	Description
115 g N	Normal net value as result
MODE: CHECKWEIGH	Activated application
UNDER LIMIT -20 pcs	Lower limit
OVER LIMIT -60 pcs	Upper limit

Check Counting application printout (Counting initialized)

Printout		Description
115 pcs	OVER	Positive net number > upper limit
- 115 pcs	OVER	Negative net number < negative upper limit
30 pcs	ACCEPT	Positive net number in the target range
15 pcs	UNDER	Positive net number < lower limit
MODE: CHECKCOU	JNT	Activated application
UNDER LIMIT -	-20 pcs	Lower limit
OVER LIMIT -	-60 pcs	Upper limit
WREF 0.35423	ΟZ	Sample weight, as displayed

Totalizing application printout (= weight total printout)

Printout		Description
N: 2		Number of items, here: 2
TOTAL:	1.955 g	Totalized value
MAX:	1.485 g	Maximum value
MIN:	0.470 g	Minimum value

Printout		Description
N: 25		Number of items, here: 25
TOTAL:	248 g	Totalized value of the weights
	124 pcs	Totalized value of the parts
MAX:	22 g	Maximum value as weight
	11 pcs	Maximum value as quantity
MIN:	4 g	Minimum value as weight
	2 pcs	Minimum value as quantity

Totalizing application printout (= parts and weight total printout)

10.2 FCC notice

Note:

This device has been tested and found to comply with the limits for digital devices of class B as per part 15 of the FCC regulations. These limits were created in order to ensure appropriate protection against interference when operating in residential areas. This device generates, uses, and may emit high-frequency energy and, if it is not installed and used in accordance with the operating instructions, may cause interference with radio communication. However, there is no guarantee that interference will not occur in certain facilities. If this device causes interference with the radio or television reception, which can be determined by switching the device off and then back on again, we recommend one or more of the following measures to eliminate the interference:

- Realignment or repositioning of the reception antenna
- Increasing the distance between the device and the receiver
- Connecting the device and the receiver to separate electric circuits
- Call in the dealer or an experienced radio/television technician

Published by Minebea Intec Bovenden GmbH & Co. KG | Leinetal 2 | 37120 Bovenden, Germany Phone: +49.551.309.83.0 | Email: info@minebea-intec.com www.minebea-intec.com

